

ELECTROLAB

Your Quality, Our Assurance

EDI-2




Bathless Disintegration Testers

Bathless Disintegration Testers

2 Station Unit



Features:

- Bathless heating saves:
 - Time
 - Electricity 
 - Water
- Quick heating of media attains media temperature 50% faster as compared to disintegration tester with Bath
- Bathless disintegration tester enhances visual observation
- Parking of baskets at the top position at the end of the test with an alarm
- Two individual sensors measure beaker temperature
- Green illumination from bottom for better visibility of disintegration process
- Programmable timer with dual range in min-sec, hrs-min and in continuous mode
- Snap-click™ loading and a swivel free basket movement
- Maintenance free smooth drive
- Suitable for veterinary tablets and support for bolus basket



SS Body

Single Time
Single Motor



654 Petrolia Road,
Toronto, ON M3J 2W3

Tel: (416) 736-6166
Toll-Free: (800) 387 3570
Fax: (416) 736-9346

Bathless Disintegration Testers

Standard Accessories



Guided disc for six tube basket assembly (Set of 6)
Part no - 0201A00004



Six tube basket assembly with guided disc (Set of 2) for EDI-2
Part no - 0202A00009



Glass beaker (1000 mL)
Part no - 0213A00001



External temperature probe for EDI-2
Part no - 0213B00001

Optional Accessories



Magnetic guided disc (Set of 12)
Part no - 0201A00007



Three tube bolus basket assembly with guided disc (Set of 2) for EDI-2
Part no - 0202A00010



Disc handling device
Part no - 0201A00018



Printer Epson LX-310 plus
Part no - 0201A00020

Bathless Disintegration Testers

Specifications

Models	EDI-2
No. of stations	2
No. of strokes	30 ± 0.5 strokes/min
Stroke height	55 mm ± 2 mm
Controlled temperature range	20°C to 39.9°C
Temperature resolution	± 0.1°C
Timer	Range 1: 1 sec to 99 min and 59 sec Range 2: 1 min to 9 hrs 59 min
Display	4 digit LED
Printer interface	—
PC connectivity	—
No. of motors	1
Illumination	Green Light
Dimensions (L x W x H)	400 mm x 365 mm x 510 mm
Weight (Kg)	18

Our Products

- Complete range of Dissolution Testers
- Dissolution Media Preparator
- Friability Testers
- Tablet Hardness Testers
- Electromagnetic Sieve Shakers
- Tap Density Tester
- Bulk Density Tester
- Powder Flow Tester
- Leak Testers
- Peristaltic Pumps

The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice
November, 2017



654 Petrolia Road,
Toronto, ON M3J 2W3

Tel: (416) 736-6166
Toll-Free: (800) 387 3570
Fax: (416) 736-9346