

ELECTROLAB

Bulk Density Tester




BETATEK INC.
Instrumentation Superbly Supported
L'instrumentation à son meilleur

654 Petrolia Road,
Toronto, ON M3J 2W3
Tel: (416) 736-6166
Toll-Free: (800) 387 3570
Fax: (416) 736-9346

Tap Density Tester




Optional Accessories



Acoustic cabinet

This reduces the sound level to 71 dB to meet laboratory standards



Adapter for 50 ml cylinder



Adapter for 10 ml cylinder

Bulk Density Tester (EV-02)



Compliance



Spatula removes excess powder

Key Features

- Complies with USP, Ph. Eur. and ASTM specifications
- Measures the bulk density of fine powders and similar products
- Easy calculation of bulk density in terms of grams per ml
- Useful in calculating powder flow-ability index