

Graunule Friability Tester



654 Petrolia Road, Toronto, ON M3J 2W3 Tel: (416) 736-6166 Toll-Free: (800) 387 3570

Fax: (416) 736-9346

Friability Testers

Dual Drum Friability Testers

- Supports 2 drums with unique interlocking mechanism to ensure correct drum position
- Patented automatic sample discharge
- Loading of samples into the drum without removing the drum
- Unique front loading system of drum
- 10° tilt facility as per USP
- Supports standard and abrasion drum
- Programmable counter and timer
- Results like % friability and weight-loss can be calculated

Key Features of EF-2W

- Variable RPM (20 to 50)
- Alphanumeric keypad to enter the product and sample information
- Stores test protocols for up to 10 different drugs
- User security with password
- Interfaced to Mettler, Sartorius and Contech weighing scales
- Serial port for complete test report printout

Granule Friability Tester (EGF-1)





Shaking

Angle

Key Features

- Complies to Ph. Eur. specifications
- Quantifiable friability of hard tablets, granules and pellets
- Bi-directional horizontal shaking action and oscillation at 42° as per USP

C€

- Contains of 105 ml glass container
- Clear acrylic lid with magnetic interlock for safe operation
- Interchangeable glass sample containers for rapid throughput
- Oscillation frequency verification certificate (optional)

